

Vacuum Controller IG-100VC



IGENE LABSERVE

Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

FEATURES

- Start and leave with evaporation detect and auto adjustments.
- GUI control panel, what you see is what you get.
- Leave your glove on, Pilot can still read your taps.
- Can be installed in every vacuuming sources.
- Suitable for laboratory or industrial environments.

SPECIFICATIONS

Model	IG-100VC
Power input	DC24V, 1A
Resolution	1 digit
Pressure display	mbar, hPa, torr
Display	7" touch screen
Measuring range	1100 ~ 0 mbar abs
Measurement accuracy	± 2 mbar (after adjustment, constant temp.)
Setting range	1013 ~ 0 mbar
Housing material	Plastic
Venting valve	Yes
Auto-cleaning	Yes
Power	15 W
Hose barb	ID 8 ~ ID10
Safety	Over-pressure Protection, Error Message
Dimension (L x W x H)	21 x 20 x 18 cm
Net weight	2.2 Kg

