

# White Light Transilluminator

**IG-WLT**



**IGENE LABSERVE**

Innovative • Interactive • Intuitive

# FEATURES & SPECIFICATIONS

## DESCRIPTION

The iGene white light transilluminator is a very cost-effective product for viewing colorimetric/protein gels and other colorimetric samples that require a white light transilluminator. It is very compact, convenient and safe to use. Its light source is LED and the lifetime bulb expectancy is 50,000 hours and power consumption are very low. It carries a one-year warranty.

## FEATURES

- ❖ **Model No.** : WLT
- ❖ **Viewing area** : 20x20cm
- ❖ **Tube** : 6
- ❖ **Light Source**: LED
- ❖ **Housing** : Anti Corrosion treated with Powder coated Mild Steel
- ❖ **Safety Shield** : Clear Transparent lid for excellent clarity
- ❖ **Input** : 8 Watt
- ❖ **Weight (kg)**: 1.3
- ❖ **Warranty** 1 year

