

Power Supply

IG-403PS



IGENE LABSERVE

Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

DESCRIPTION

Power supply comes with a LCD screen display which displays the voltage, current, power, and timing information and has a programmable storage up to 10 programs

FEATURES

- ✦ LCD screen displays voltage, current, programs, Pause and timing information
- ✦ Output voltage ranges 0-300V dc
- ✦ Current ranges 1mA to 999mA
- ✦ It comes with 10 programmable memory storage
- ✦ It has Pause key to pause the running program
- ✦ It recalls the parameter after power failure



SPECIFICATIONS

Model No.	IG-403PS
Output voltage	0-300V, fully adjustable in 1 V step
Current	1mA to 999 mA
Watt	300W
Digital	Digital LCD display
Output terminals	4 pair recessed banana jacks in parallel
Time	1 min to 99 mins
Programs	10 program memory
Operating conditions	0-40°C; 0-95% humidity
Safety features	No- load detection, overload/short circuit detection, overvoltage protection, over temp protection

